

CBS™ High Security Straw

The only secure system
for LN₂ storage



Sales and Marketing

29, rue Tronchet
75008 PARIS - FRANCE
Tel. +33 (0) 1 49 24 05 05
Fax +33 (0) 1 49 24 05 01

e-mail : bl@cryobiosystem-imv.com



www.cryobiosystem.com

Order Processing

10, rue Clemenceau
61300 L'AIGLE - FRANCE
Tel. +33 (0) 2 33 34 64 64
Fax +33 (0) 2 33 34 11 98

e-mail : cv@cryobiosystem-imv.com

CE Mark for Medical Device CE 93/42 F.D.A. 510 (k) clearance



CE 93/42
EN 46002
ISO 9002
UL Certification
FDA 510 (k) clearance

- Sterilized by gamma irradiation
- Blister compartments sealed with Tyvek®
- ID labels and jackets (with special LN₂ resistant labels):
 - Easy to use
 - Clearly identified text
 - Flexibility of text to be written
 - Flexible barcoded and/or alphanumerical text
- ID rod offers color code (6 colors are available)
- Tamperproof identification:
 - Rod is inserted inside the sealed CBS™ straw,
 - With a hydrophobic plug safely separating the embryo from the ID rod
- Hydrophobic plug for embryo straw:
 - New patented hydrophobic, non absorbing plug
 - No absorption of media or cells inside the plug
- Biocompatibility tests:
 - Straw/content physiochemical inertia
 - = non toxic (see technical product sheet)
- Security: Heat sealed CBS™ High Security straw is absolutely leak-proof in LN₂ and mechanically resistant to pressures of up to 150 kg/cm²
- Increase storage capacity
(3 to 5 times more storage than with vials)

Product Code	Product	Packaging / Sales Unit
010286	CBS™ Embryo Straw	1 straw/compartament 5 blisters x 4 compartments
010266 to 271	ID Rod (available in 6 colors)	In bags of 10 per color
007214	SYMS Pump	Per unit
016399	SYMS Sealing unit 220 Volt - 50 Hz	Per unit
017060	SYMS Sealing unit 110 Volt - 60 Hz	Per unit
010287	CBS™ Sperm Straw	5 straws/compartament
010288	CBS™ Sperm Straw	20 straws/compartament